

POWER MOBILITY CHART



	Cobalt Standard/Transportable	Titan Standard/Transportable	Titan P22 Standard	Image EC Standard	Trident Standard	Trident HD Heavy-Duty	Cirrus Plus EC Folding	Cirrus Plus HD Heavy-Duty Folding
HCPCS	K0821/K0823	K0816 / K0821 / K0823	K0822	K0823	K0823	K0825	K0822	K0824
Weight Capacity	250 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	450 lbs.	300 lbs.	400 lbs.
Top Speed	4 mph	4 mph	4 mph	4 mph	4 mph	4 mph	5 mph	5 mph
Maximum Range	8 miles	15 miles	15 miles	15 miles	15 miles	15 miles	15 miles	15 miles
Turning Radius	29"	25"	25"	27.5"	24"	29"	30" 32" 33"	35"
Climbing Angle	6°	6°	6°	6°	7.5°	6°	6°	6°
Ground Clearance	2"	2.5"	2.5"	2"	2"	2.5"	4"	4"
Dimensions (L x W)	38.5" x 24"	40" x 24"	40" x 24"	42" x 24"	42.25" x 27.25"	48" x 29"	42.5" x 24" 42.5" x 26" 42.5" x 28"	42.5" x 29"
Floor-to-Seat Height	20" - 23"	20.5" - 23"	20.5" - 23"	20" - 23"	21" - 24"	21" - 23"	18.5"	20"
Seat Type	Stadium	Captain's	Pan Seat	Captain's	Captain's	Captain's	Sling	Sling
Seat Dims. (WxD)	18" x 15"	16"x16" 16"x18" 18"x16" 18"x18" 20"x18" 20"x20" 22"x20"	18" x 18" / 20" x 18"	18" x 18" 20" x 18"	18" x 18" 20" x 18"	22" x 20" 24" x 20"	16" x 16" 18" x 16" 20" x 16"	22" x 16"
Controller	Dynamic Shark 40A	Dynamic Shark 40A	Dynamic Shark 40A	PG VR2 50A	PG VR2 50A	PG VR2 70A	PG VSi 50A	PG VSi 50A
Motor	24V x 120W x 4700 rpm	24V x 200W x 3900 rpm	24V x 200W x 3900 rpm	24V x 200W x 3900 rpm	24V x 320W x 4600 rpm	24V x 450W x 2700 rpm	24V x 320W x 4600 rpm	24V x 320W x 4600 rpm
Batteries (pr)	12V x 21AH	12V x 35AH (U1)	12V x 35AH (U1)	12V x 35AH (U1)	12V x 36AH (U1)	12V x 50AH	12V x 36AH (U1)	12V x 36AH (U1)
Charger	2A Offboard	3A Offboard	3A Offboard	4A Offboard	4A Offboard	5A Offboard	4A Offboard	4A Offboard
Brakes	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic
Front Wheels	8" x 2" Solid	6" x 2" Solid	6" x 2"	6" x 2" Solid	8" x 2" Solid	10" x 3" Solid	8" x 2" Solid	8" x 2" Solid
Rear Wheels	9" x 3" Flat-free	10" x 3" Flat-free	10" x 3" Flat-free	10" x 3" Flat-free	14" x 3" Flat-free	14" x 3" Flat-free	12.5" x 2.25" Flat-free	12.5" x 2.25" Flat-free
Color	Red or Blue	Red & Blue	Red & Blue	Red & Blue	Black	Black	Black	Black
Base Weight	74 lbs.	68 lbs.	68 lbs.	88 lbs.	82 lbs.	98 lbs.	90 lbs.	90 lbs.
Battery Weight (pr)	28 lbs.	44 lbs.	44 lbs.	53 lbs.	53 lbs.	66 lbs.	56 lbs.	56 lbs.
Seat Weight	19 lbs.	43 lbs.	43 lbs.	51 lbs.	51 lbs.	48 lbs.	N/A	N/A
Total Weight	121 lbs.	158 lbs.	158 lbs.	192 lbs.	186 lbs.	212 lbs.	146 lbs.	146 lbs.
Warranty on Frame	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime
Warranty on Electronic/ Drive Train Components	14 months	14 months	14 Months	14 Months	14 Months	14 Months	14 Months	14 Months
Warranty on Batteries	6 months	6 months	6 Months	6 Months	6 Months	6 Months	6 months	6 months
MSRP	\$5,295.00	\$5,560.00 (18" SB) \$5,725.50 (18" HB) \$5,891.00 (20")	\$5,560.00 (18") \$5,725.50 (20")	\$6,414.00 (18") \$6,612.00 (20")	\$6,387.00 (18") \$6,586.00 (20")	\$9,094.00 (22") \$9,260.00 (24")	\$6,222.00 (16"/18") \$6,784.50 (20")	\$8,439.00