

# Titan X23 Front Wheel Drive Power Wheelchair Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>	<b>P.O. #</b>	<b>Date:</b>	
<b>Bill to:</b>		<b>Ship to:</b>	

## Powerchair with Seating

### Standard Back Captain's Seat

- TITANLB18CSX23 .....18" (W) x 18" (D) .....\$5,560.00  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN1616X23 .....16" (W) x 16" (D) .....\$5,395.00  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN1618X23 .....16" (W) x 18" (D) .....\$5,460.50  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN1816X23 .....18" (W) x 16" (D) .....\$5,494.00  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*

### High Back Captain's Seat

- TITAN18CSX23 .....18" (W) x 18" (D) .....\$5,725.50  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN20CSX23 .....20" (W) x 18" (D) .....\$5,891.50  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN2020X23 .....20" (W) x 20" (D) .....\$5,957.00  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*
- TITAN22CSX23 .....22" (W) x 20" (D) .....\$6,056.00  
*HCPCS: K0823, Includes interchangeable color shrouds in Blue & Red.*

## Standard Armrests

- ARM-L.....Included  
*HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests*
- ARM-R .....Included  
*HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests*

## Controller Arm

- Adjustable Length.....Included
- AA4100 – Swing-Away Controller Arm (HCPCS E1028) .....\$105.00

## Front Riggings

- Adjustable Footplate .....Included
- AE2500/LK3JELR – Elevating Legrest Kit .....\$207.48  
*HCPCS: E0990, Includes swing-away elevating legrest with attachment bracket*

## Batteries & Charger

- U1BATTERY (pair) .....Included  
*HCPCS: E2365, 12V x 35AH*
- 3A Battery Charger .....Included

## Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978) .....Included
- AH1000 Crutch/Cane Holder (HCPCS E2207) .....\$120.12
- AH2000 Cylinder Tank Carrier (HCPCS E2208) .....\$156.52
- ST205-COVER 18" Seat Cover .....\$118.17
- ST301-COVER 20" Seat Cover .....\$118.17
- AZ0060 Cup Holder .....\$50.96

## Specifications

<b>HCPCS</b>	
TITANLB18CSX23	K0823
TITAN1616X23	K0823
TITAN1618X23	K0823
TITAN1816X23	K0823
TITAN18CSX23	K0823
TITAN20CSX23	K0823
TITAN2020X23	K0823
TITAN22CSX23	K0823
<b>Performance</b>	
Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	25"
<b>Dimensions</b>	
Overall Length	40"
Overall Width	25"
Floor to Seat Height	20.5"–23"
Ground Clearance	2.5"
<b>Power</b>	
Motor	Two 24V x 200W x 3900rpm
Controller	Dynamic Shark 40A
Brakes	Electromagnetic
Freewheel Mode	Yes
<b>Tires</b>	
Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free
<b>Weights</b>	
Base Weight	68 lbs.
Battery Weight	44 lbs.
Seat Weight	43 lbs.
Heaviest Piece Weight (Disassembled)	54 lbs.
Total Weight	155 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.