

**Trident HD Front Wheel Drive Heavy Duty
Power Wheelchair
Order Form/Price List**



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Powerchair with Seating

High Back Captain's Seat

- 2850HD-2222" (W) x 20" (D)\$5,001.36
HCPCS: K0825
- 2850HD-2424" (W) x 22" (D)\$5,092.36
HCPCS: K0825

Standard Armrests

- ARM-L.....Included
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests
- ARM-RIncluded
HCPCS: E0973, Adjustable Height, Removable, Flip-Up Armrests

Controller Arm

- Adjustable Length.....Included
- TRID-2A37 – Swing-Away Controller Arm (HCPCS E1028)\$105.00

Front Riggings

- Adjustable FootplateIncluded
- TRID-2A47/SFELR/LK3JELR – Elevating Legrest Kit.....\$207.48
HCPCS: E0990, Includes swing-away elevating legrest with attachment bracket

Batteries & Charger

- 22NFBATTERY (pair)Included
HCPCS: E2361, 12V x 50AH
- 5A Battery ChargerIncluded

Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978)Included
- AH1000 Crutch/Cane Holder (HCPCS E2207)\$120.12
- AH2000 Cylinder Tank Carrier (HCPCS E2208).....\$156.52
- AZ0060 Cup Holder\$50.96

Specifications

HCPCS

2850HD-22	K0825
2850HD-24	K0825

Performance

Maximum Weight Capacity	450 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	29"

Dimensions

Overall Length	48"
Overall Width	29.5"
Floor to Seat Height	21"–23"
Ground Clearance	2.5"

Power

Motor	Two 24V x 450W x 2700rpm
Controller	PG VR2 70A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires

Casters	10" x 3" Solid
Drive Wheels	14" x 3" Flat Free

Weights

Base Weight	98 lbs.
Battery Weight	66 lbs.
Seat Weight	48 lbs.
Total Weight	212 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.