

# Cobalt X23 Travel Rear Wheel Drive Power Wheelchair

## Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>	<b>P.O. #</b>	<b>Date:</b>	
<b>Bill to:</b>		<b>Ship to:</b>	

### Chair with Seating

- Stadium Seat**
- COBALTX23BL16FS ..... Midnight Blue.....\$5295.00
  - COBALTX23RD16FS ..... Metallic Red .....\$5295.00

### Armrests

- Removable, Full Length, Flip-up Armrests (HCPCS E0973) ..... Included

### Controller Arm

- Adjustable Length ..... Included
- AA4100 Swingaway Controller Arm (HCPCS E1028) .....\$105.00

### Front Riggings

- Adjustable Footplate ..... Included
- AE2500/LK3JELR Elevating Legrest Kit .....\$328.00  
Includes left & right Swingaway ELRs (HCPCS E0990) & attachment bracket.

### Batteries & Charger

- 12V x 21AH Batteries (pair) (HCPCS K0733) ..... Included
- 2A Battery Charger ..... Included

### Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978) ..... Included
- AH1000 Crutch/Cane Holder (HCPCS E2207) .....\$53.00
- AH2000 Cylinder Tank Carrier (HCPCS E2208).....\$70.00
- AB1000 Armrest Bag .....\$25.00
- AB1120 Backpack .....\$56.00
- AZ0060 Cup Holder .....\$30.00

### Specifications

<b>HCPCS</b>	<b>K0823</b>
<b>Maximum Weight Capacity</b>	250 lbs.
<b>Maximum Speed</b>	4 mph
<b>Maximum Range</b>	7 miles
<b>Maximum Grade/Incline</b>	6°
<b>Turning Radius</b>	29"
<b>Ground Clearance</b>	2"
<b>Dimensions</b>	
Overall Length	38.5"
Overall Width	24"
Floor to Seat Height	20"-23"
<b>Power</b>	
Motor	Two 24V x 120W x 4700rpm
Controller	Dynamic A50A
Batteries	Two DC12V x 21AH
Battery Charger	2 Amp, Offboard
<b>Brakes</b>	Electromagnetic
<b>Freewheel Mode</b>	Yes
<b>Tires</b>	
Front Casters (Solid)	8" x 2"
Rear Drive Wheels (Flat Free)	9" x 3"
<b>Weights</b>	
Without Batteries	93 lbs.
With Batteries	121 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.*