

**Image EC Mid Wheel Drive Standard  
Power Wheelchair  
Order Form/Price List**



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>		<b>P.O. #</b>	<b>Date:</b>
<b>Bill to:</b>		<b>Ship to:</b>	

**Chair with Seating**

**Captain's Seat**

- 2800ECBU-RCL..... 18" (W) x 18" (D).....\$6,414.00  
*Includes interchangeable color shrouds in Blue & Red.*
- 2800ECBU-RCL-20 ..... 20" (W) x 18" (D) .....\$6,612.00  
*Includes interchangeable color shrouds in Blue & Red.*

**Armrests**

- Adjustable Height, Removable, Full Length, Flip-up Armrests (HCPCS E0973).....Included

**Controller Arm**

- Adjustable Length..... Included
- SF1009 Swingaway Controller Arm (HCPCS E1028).....\$247.50

**Front Riggings**

- Adjustable Footplate ..... Included
- SFELR/LK3JELR/ECELR-P Elevating Legrest Kit .....\$328.00  
*Includes left & right Swingaway ELRs (HCPCS E0990) & attachment bracket.*

**Batteries & Charger**

- 12V x 35AH (U1) Batteries (pair) (HCPCS E2365) ..... Included
- 4A Battery Charger ..... Included

**Accessories**

- Positioning/Safety Belt (HCPCS E0978) ..... Included
- AH1000 Crutch/Cane Holder (HCPCS E2207) .....\$53.00
- AH2000 Cylinder Tank Carrier (HCPCS E2208).....\$70.00
- AB1000 Armrest Bag .....\$25.00
- AB1120 Backpack .....\$56.00
- ST205-COVER 18" Seat Cover .....\$118.17
- ST301-COVER 20" Seat Cover .....\$118.17
- AZ0060 Cup Holder .....\$30.00

**Specifications**

<b>HCPCS</b>	K0823
<b>Performance</b>	
Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	27.5"
<b>Dimensions</b>	
Overall Length	42"
Overall Width	24"
Floor to Seat Height	20"-23"
Ground Clearance	2"
<b>Power</b>	
Motor	Two 24V x 200W x 3900rpm
Controller	PG VR2 50A
Brakes	Electromagnetic
Freewheel Mode	Yes
<b>Tires</b>	
Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free
<b>Weights</b>	
Base Weight	88 lbs.
Battery Weight	53 lbs.
Seat Weight	51 lbs.
Total Weight	192 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.