

# Spitfire Scout 4-Wheel Transport Scooter

## Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>		<b>P.O. #</b>	<b>Date:</b>
<b>Bill to:</b>		<b>Ship to:</b>	

### Scooter with Seating

#### Folding Seat

- SFSCOUT4.....16.5" (W) x 16" (D) .....\$1,820.00  
*Includes interchangeable color shrouds in Blue & Red.*
- SFSCOUTDLX4.....Deluxe Model, 16.5" (W) x 16" (D) .....\$1,874.50  
*Includes interchangeable color shrouds in Blue & Red.*

### Batteries & Charger

- 12V x 12AH Batteries (pair) ..... Included
- 2A Battery Charger ..... Included

### Accessories

- AB1000 Armrest Bag .....\$25.00
- AB1120 Backpack .....\$56.00
  
- AH1000 Crutch/Cane Holder .....\$53.00
- AH2000 Cylinder Tank Carrier .....\$70.00
- AB2100 Rear Basket .....\$70.00
  
- AZ1000 Scooter Cover .....\$62.00
  
- AT1000 Scooter Trailer .....\$204.00
  
- AZ0060 Cup Holder .....\$30.00

### Specifications

**HCPCS** **K0800 (SFSCOUT4 only)**

#### Performance

Maximum Weight Capacity	300 lbs.
Maximum Speed	4.25 mph
Maximum Range	9 miles
Maximum Grade/Incline	6°
Turning Radius	53.75"

#### Dimensions

Overall Length	42.5"
Overall Width	21.25"
Floor to Seat Height	17"-19"
Ground Clearance	2.5"

#### Power

Motor	24V x 270W x 4700rpm
Controller	CTE 45A
Brakes	Electromagnetic
Freewheel Mode	Yes

#### Tires (Flat-free)

Front	8" x 2"
Rear	8" x 2"

#### Weights

Base Weight	58 lbs.
Battery Weight	20 lbs.
Seat Weight	16 lbs.
Total Weight	94 lbs.
Heaviest Piece Weight	39 lbs. (Front Section)

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*