

Titan Front Wheel Drive Standard/Transportable Power Wheelchair Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Chair with Seating

Standard Captain's Seat

- TITANLB18CS18" (W) x 18" (D)\$5,560.00
Includes interchangeable color shrouds in Blue & Red.
- TITANLB18CSX2318" (W) x 18" (D)\$5,560.00
Includes interchangeable color shrouds in Blue & Red.

High Back Captain's Seat

- TITAN18CS18" (W) x 18" (D)\$5,725.50
Includes interchangeable color shrouds in Blue & Red.
- TITAN18CSX2318" (W) x 18" (D)\$5,725.50
Includes interchangeable color shrouds in Blue & Red.
- TITAN20CS20" (W) x 18" (D)\$5,891.00
Includes interchangeable color shrouds in Blue & Red.
- TITAN20CSX2320" (W) x 18" (D)\$5,891.00
Includes interchangeable color shrouds in Blue & Red.

Armrests

- Adjustable Height, Removable, Full Length, Flip-up Armrests
(HCPCS E0973) Included

Controller Arm

- Adjustable Length Included
- AA4100 Swingaway Controller Arm (HCPCS E1028)\$105.00

Front Riggings

- Adjustable Footplate Included
- AE2500/LK3JELR Elevating Legrest Kit\$207.48
Includes left & right Swingaway ELRs (HCPCS E0990) & attachment bracket.

Batteries & Charger

- 12V x 35AH (U1) Batteries (pair) (HCPCS E2365) Included
- 3A Battery Charger Included

Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978) Included
- AH1000 Crutch/Cane Holder (HCPCS E2207)\$120.12
- AH2000 Cylinder Tank Carrier (HCPCS E2208)\$156.52
- AB1000 Armrest Bag\$52.78
- AB1120 Backpack\$163.80
- ST205-COVER 18" Seat Cover\$118.17
- ST301-COVER 20" Seat Cover\$118.17
- AZ0060 Cup Holder\$50.96

Specifications

HCPCS	
TITANLB18CS/TITAN18CS/TITAN20CS	K0821
TITANLB18CSX23/TITAN18CSX23/TITAN20CSX23	K0823
Performance	
Maximum Weight Capacity	300 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	25"
Dimensions	
Overall Length	40"
Overall Width	25"
Floor to Seat Height	20.5"-23"
Ground Clearance	2.5"
Power	
Motor	Two 24V x 200W x 3900rpm
Controller	Dynamic Shark 40A
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires	
Casters	6" x 2" Solid
Drive Wheels	10" x 3" Flat Free
Weights	
Base Weight	68 lbs.
Battery Weight	44 lbs.
Seat Weight	43 lbs.
Total Weight	155 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.