

Phoenix HD 3-Wheel Heavy Duty Transport Scooter

Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Stadium Seat

- PHOENIXHD3 17.5" (W) x 17" (D)\$2,548.00
Includes interchangeable color shrouds in Blue & Red.

Batteries & Charger

- 12V x 20AH Batteries (pair) Included
- 1.5A Battery Charger Included

Accessories

- AB1000 Armrest Bag\$25.00
- AB1120 Backpack\$56.00
- AH1000 Crutch/Cane Holder\$53.00
- AH2000 Cylinder Tank Carrier\$70.00
- AB2100 Rear Basket\$70.00
- AZ1000 Scooter Cover\$62.00
- AT1000 Scooter Trailer\$204.00
- AZ0060 Cup Holder\$30.00

Specifications

HCPCS **K0801**

Performance

Maximum Weight Capacity	350 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	43"

Dimensions

Overall Length	39"
Overall Width	22.5"
Floor to Seat Height	20"-22"
Ground Clearance	3.75"

Power

Motor	24V x 350W x 3600rpm
Controller	CTE 70A
Brakes	Electromagnetic
Freewheel Mode	Yes

Tires (Flat-free)

Front	8" x 2"
Rear	8" x 2"

Weights

Base Weight	63 lbs.
Battery Weight	31 lbs.
Seat Weight	21 lbs.
Total Weight	115 lbs.
Heaviest Piece Weight	32 lbs. (Front Section)

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.