

SCOOTER MOBILITY CHART



	Bobcat X Transport	ZooMe Flex Transport	Scout Transport	Scout DST Transport	Phoenix HD Heavy-Duty Transport	Ventura Mid-Size	Ventura DLX Mid-Size	Prowler Executive Scooter	Odyssey GT Executive Scooter	Cobra GT4 Executive Scooter	King Cobra Executive Scooter	King Cobra PGV Executive Scooter	Maverick Executive Scooter
HPCPS	N/A	N/A	K0800	N/A	K0801	K0807	N/A	K0899	N/A	N/A	N/A	N/A	N/A
Weight Capacity	265 lbs.	275 lbs.	300 lbs.	325 lbs.	350 lbs.	350 lbs.	350 lbs.	500 lbs.	325 lbs.	450 lbs.	450 lbs.	450 lbs.	400 lbs.
Top Speed	4 mph	4 mph	4.25 mph	4.25 mph	4 mph	6 mph	6 mph	7 mph	9.5 mph	10 mph	11 mph	10 mph	9.5 mph
Maximum Range	7.5 miles	13 miles	9 miles / 15 miles	15 miles	15 miles	22 miles	22 miles	23 miles	22 miles	35 miles	35 miles	35 miles	35 miles
Turning Radius	45" (3) 53" (4)	53"	45.5" (3) 53.75" (4)	45" (3) / 55" (4)	43" (3) 54" (4)	53" (3) 54" (4)	53" (3) 54" (4)	50" (3) 69" (4)	47.25"	82"	83"	83"	62"
Climbing Angle	6°	6°	6°	6°	6°	7.5°	7.5°	7.5°	6°	8°	8°	8°	10°
Ground Clearance	2.5"	2"	2.5"	2.5"	3.75"	2.5"	2.5"	5"	3.95"	7"	4"	4"	5"
Dimensions (L x W)	42" x 19.25"	37" x 18"	42.5" x 20.5"	42.5"x 21.25"	41.5" x 22.5"	48" x 25"	48" x 25"	51" x 26"	52"x 26.5"	63" x 29"	67" x 28.8"	85" x 29"	62" x 29"
Floor-to-Seat Height	17" - 19"	19.5"	17" - 19"	17" - 19"	20" - 22"	21.5" - 23.5"	21.5" - 23.5"	25" - 28"	21.5" - 23.5"	27" - 30"	27" - 30"	27" - 30"	27"
Seat Type	Stadium	Stadium	Stadium	Stadium	Stadium	Stadium	Captain's	Captain's	Sport Style Captain's	Sport Style Captain's	Sport Style Captain's	Sport Style Captain's	Sport Style Captain's
Seat Dims. (WxD)	16.5" x 13.5"	16.75" x 15"	16.5" x 16"	17.5" x 16"	17.5" x 17", 20" x 17"	18" x 17.5"	18" x 18" 20" x 18"	20" x 20" 22" x 20"	18" x 18"	20" x 20" 22" x 20"	22" x 20"	22" x 20"	20" x 20" 22" x 20"
Controller	PG S Drive 45A	PG S Drive 45A	Curtis 45A	Curtis 45A	CTE 70A	PG S Drive 70A	PG S Drive 70A	Dynamic Rhino II 120A	PG S DRIVE 120A	Dynamic Rhino II 160A	Dynamic Rhino II 180A	Dynamic Rhino II 180A	Dynamic Rhino II 160A
Motor	24V x 270W x 4700 rpm	24V x 120W	24V x 270W x 4700 rpm	24V x 270W x 4700 rpm	24V x 350W x 3600 rpm	24V x 350W x 5100 rpm	24V x 350W x 5100 rpm	24V x 1300W x 4350 rpm	24V x 500W x 5100 rpm	24V x 950W x 5200 rpm	24V x 1300W x 5200 rpm	24V x 1300W x 5200 rpm	24V x 1300W x 4300 rpm
Batteries (pr)	12Vx10AH	25.9Vx10AH Lithium	12Vx12AH / 12Vx21AH	12Vx21AH	12Vx20AH	12Vx33AH	12Vx33AH	12V x 50AH	12V x 50AH	12V x 75AH	12V x 75AH	12V x 75AH	12V x 75AH
Charger	1.5A Offboard	2A Offboard	2A Offboard	2A Offboard	1.5A Offboard	4A Offboard	4A Offboard	5A Offboard	5A Offboard	5A Offboard	8A Offboard	8A Offboard	8A Offboard
Brakes	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic	Electromagnetic & Hand	Electromagnetic	Electromagnetic & Hand	Electromagnetic & Hand	Electromagnetic & Hand	Electromagnetic & Hand
Front Wheels	8" x 2" Flat-Free	6" x 2"	8" x 2" Flat-Free	8" x 2"	9" x 3"	10" x 3" Flat-Free	10" x 3" Flat-Free	13" x 5" Pneumatic	11" x 3.5" Pneumatic	14" x 4" Pneumatic	16.25" x 4" Pneumatic	14.5" x 6"	19.5" x 2.75" Pneumatic
Rear Wheels	8" x 2" Flat-Free	7" x 2"	8" x 2" Flat-Free	8" x 2"	9" x 3"	10" x 3" Flat-Free	10" x 3" Flat-Free	13" x 5" Pneumatic	11" x 4" Pneumatic	14" x 4" Pneumatic	17" x 4" Pneumatic	17" x 6"	19.5" x 2.75" Pneumatic
Color	Red	Blue	Red & Blue	Red & Blue	Red & Blue	Red & Blue	Red & Blue	Metallic Grey	Electric Blue	Metallic Grey	Vivid Red	Vivid Red	Black
Base Weight	50 lbs. (3) 59 lbs.(4)	53.6 lbs.	49 lbs. (3) 58 lbs. (4)	60 lbs. (3) / 70 lbs. (4)	63 lbs. (3) 70 lbs. (4)	88 lbs. (3) 99 lbs. (4)	88 lbs. (3) 99 lbs. (4)	130 lbs. (3) 140 lbs. (4)	133 lbs.	191 lbs.	209 lbs.	209 lbs.	210 lbs.
Battery Weight (pr)	18 lbs.	5 lbs.	20 lbs.	30 lbs.	31 lbs.	44 lbs.	44 lbs.	66 lbs.	66 lbs.	94 lbs.	94 lbs.	94 lbs.	106 lbs.
Seat Weight	16 lbs.	N/A	16 lbs.	23 lbs.	21 lbs.	26 lbs.	55 lbs.	55 lbs.	55 lbs.	55 lbs.	60 lbs.	60 lbs.	50 lbs. 55 lbs.
Total Weight	84 lbs. (3) 93 lbs. (4)	58.6 lbs.	85 lbs. (3) 94 lbs. (4)	113 lbs. (3) / 123 lbs. (4)	115 lbs. (3) 129 lbs. (4)	158 lbs. (3) / 169 lbs. (4)	187 lbs. (3) / 198 lbs. (4)	251 lbs. (3) 261 lbs. (4)	254 lbs.	340 lbs.	363 lbs.	363 lbs.	366 lbs. 371 lbs.
Heaviest Piece Weight	28 lbs. (3) (Front Sec.) 38 lbs. (4) (Front Sec.)	53.6 lbs. (Base)	29 lbs. (3) (Front Sec.) 39 lbs. (4) (Front Sec.)	34 lbs. (3) (Front Sec.) 43 lbs. (4) (Front Sec.)	32 lbs. (3) (Front Sec.) 46 lbs. (4) (Front Sec.)	60 lbs. (3) (Front Sec.) 70 lbs. (4) (Front Sec.)	60 lbs. (3) (Front Sec.) 70 lbs. (4) (Front Sec.)	N/A	N/A	N/A	N/A	N/A	N/A
Warranty on Frame	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime	Lifetime
Warranty on Electronic/ Drive Train Components	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months	14 months
Warranty on Batteries	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months	6 months
MSRP (3-Wheel)	\$1,638.00	N/A	\$1,765.50 (12AH) / \$1,874.50 (21AH)	\$2,162.00	\$2,548.00	\$3,112.00	\$3,367.00 (18") \$3,549.00 (20")	\$5,311.00 (20") \$5,456.50 (22")	N/A	N/A	N/A	N/A	\$10,556.00 (20") \$10,738.00 (22")
MSRP (4-Wheel)	\$1,729.00	\$4,273.00	\$1,820.00 (12AH) / \$1,929.00 (21AH)	\$2,235.00	\$2,730.00	\$3,330.50	\$3,549.00 (18") \$3,731.00 (20")	\$5,493.00 (20") \$5,602.00 (22")	\$7,069.00	\$8,645.00 (20") \$9,009.00 (22")	\$10,009.00	\$11,739.00	N/A