

Cobalt Travel Rear Wheel Drive Power Wheelchair Order Form/Price List



Company Name:		Phone:	
Contact Person:		P.O. #	Date:
Bill to:		Ship to:	

Chair with Seating

- Stadium Seat**
- COBALTLB16FS Midnight Blue.....\$5295.00
 - COBALTRD16FS Metallic Red\$5295.00

Armrests

- Removable, Full Length, Flip-up Armrests (HCPCS E0973) Included

Controller Arm

- Adjustable Length Included
- AA4100 Swingaway Controller Arm (HCPCS E1028)\$105.00

Front Riggings

- Adjustable Footplate Included
- AE2500/LK3JELR Elevating Legrest Kit\$328.00
Includes left & right Swingaway ELRs (HCPCS E0990) & attachment bracket.

Batteries & Charger

- 12V x 21AH Batteries (pair) (HCPCS K0733) Included
- 2A Battery Charger Included

Accessories

- LDR410017 Positioning/Safety Belt (HCPCS E0978) Included
- AH1000 Crutch/Cane Holder (HCPCS E2207)\$53.00
- AH2000 Cylinder Tank Carrier (HCPCS E2208).....\$70.00
- AB1000 Armrest Bag\$25.00
- AB1120 Backpack\$56.00
- AZ0060 Cup Holder\$30.00

Specifications

HCPCS	K0821
Maximum Weight Capacity	250 lbs.
Maximum Speed	4 mph
Maximum Range	7 miles
Maximum Grade/Incline	6°
Turning Radius	29"
Ground Clearance	2"
Dimensions	
Overall Length	38.5"
Overall Width	24"
Floor to Seat Height	20"-23"
Power	
Motor	Two 24V x 120W x 4700rpm
Controller	Dynamic A50A
Batteries	Two DC12V x 21AH
Battery Charger	2 Amp, Offboard
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires	
Front Casters (Solid)	8" x 2"
Rear Drive Wheels (Flat Free)	9" x 3"
Weights	
Without Batteries	93 lbs.
With Batteries	121 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.