

Company Name:			Phone:	
Contact Person:	P.O. #			Date:
Bill to:		Ship to:		

Cha	in with Continu
Cna	air with Seating
Sling	G Seat CPN16FBA EC, 16" (W) x 16" (D) \$6,222.00   CPN18FBA EC, 18" (W) x 16" (D) \$6,222.00   CPN20FBA EC, 20" (W) x 16" (D) \$6,784.50   CPN22FBA HD, 22" (W) x 16" (D) \$8,439.00
Arn	nrests
	Adjustable Height, Removable, Full Length, Flip Up Armrests
	(HCPCS E0973)Included
Cor	ntroller Arm
	Adjustable LengthIncluded
	AA4000 Swingaway Controller Arm (HCPCS E1028)\$247.50
Fro	nt Riggings
Ш	Adjustable FootplateIncluded
	CPELR Elevating Legrest Kit\$328.00
	Includes left & right swing-away ELRs (HCPCS E0990).
Bat	teries & Charger
	12V x 36AH (U1) Batteries (pair) (HCPCS E2365)Included
	4A Battery ChargerIncluded
Δ	
ACC	essories
	Positioning/Safety Belt (HCPCS E0978)
	AB1010 Armrest Bag\$25.00

Constitutions	
Specifications	
HCPCS	
Cirrus Plus EC	K0822
Cirrus Plus HD	K0824
Maximum Weight Capacity	
Cirrus Plus EC	300 lbs.
Cirrus Plus HD	400 lbs.
Performance	
Maximum Speed	5 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	
CPN16FBA	30"
CPN18FBA	32"
CPN20FBA	33"
CPN20FBA	35"
Dimensions	
Overall Length	42.5"
Overall Width (CPN16FBA)	24"
Overall Width (CPN18FBA)	26"
Overall Width (CPN20FBA)	28"
Overall Width (CPN20FBA)	29"
Floor to Seat Height	18.5"
Ground Clearance	4"
Power	
Motor	Two 24V x 320W x 4600rpm
Controller	PG VSi 50A
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires	
Casters	8" x 2" Solid
Drive Wheels	12.5" x 2.25" Flat Free
Weights	
Base Weight	90 lbs.
Battery Weight	56 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

146 lbs.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

Total Weight