

# Wildcat 450 Heavy Duty Folding Power Wheelchair Order Form/Price List



<b>Company Name:</b>		<b>Phone:</b>	
<b>Contact Person:</b>	<b>P.O. #</b>	<b>Date:</b>	
<b>Bill to:</b>		<b>Ship to:</b>	

### Chair with Seating

**Sling Seat**  
 WILDCAT450BK24SS..... 24" (W) x 20" (D) .....\$9,657.00

### Armrests

Adjustable Height, Removable, Full Length, Flip-up Armrests  
 (HCPCS E0973)..... Included

### Controller Arm

Adjustable Length..... Included  
 AA4000 Swingaway Controller Arm (HCPCS E1028) .....\$247.50

### Front Riggings

Adjustable Footplate ..... Included  
 AE2100-450 Elevating Legrest Kit .....\$328.00  
 Includes left & right Swingaway ELRs (HCPCS E0990) & attachment bracket.

### Batteries & Charger

12V x 50AH (22NF) Batteries (pair) (HCPCS E2361) ..... Included  
 4A Battery Charger ..... Included

### Accessories

Positioning/Safety Belt (HCPCS E0978) ..... Included  
 AB1010 Armrest Bag .....\$25.00

### Specifications

<b>HCPCS</b>	<b>K0824</b>
<b>Performance</b>	
Maximum Weight Capacity	450 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	41"
<b>Dimensions</b>	
Overall Length	48"
Overall Width	26.5"
Floor to Seat Height	20"
Ground Clearance	4"
<b>Power</b>	
Motor	Two 24V x 350W x 3800rpm
Controller	Dynamic Shark 60A
Brakes	Electromagnetic
Freewheel Mode	Yes
<b>Tires</b>	
Casters	9" x 2.5" Solid
Drive Wheels	12.5" x 2.25" Flat Free
<b>Weights</b>	
Base Weight	100 lbs.
Battery Weight	73 lbs.
Total Weight	173 lbs.

*Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.*

*Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.*