

Phoenix HD 4-Wheel Heavy Duty Transport Scooter

Order Form/Price List



Company Name:		Phone:	
Contact Person:		P.O. #	Date:
Bill to:		Ship to:	

Scooter with Seating

- Stadium Seat**
- PHOENIXHD4 17.5" (W) x 17" (D)\$2,730.00
Includes interchangeable color shrouds in Blue & Red.

Batteries & Charger

- 12V x 20AH Batteries (pair) Included
- 1.5A Battery Charger Included

Accessories

- AB1000 Armrest Bag\$25.00
- AB1120 Backpack\$56.00
- AH1000 Crutch/Cane Holder\$53.00
- AH2000 Cylinder Tank Carrier\$70.00
- AB2100 Rear Basket\$70.00
- AZ1000 Scooter Cover\$62.00
- AT1000 Scooter Trailer\$204.00
- AZ0060 Cup Holder\$30.00

Specifications

HCPCS	K0801
Performance	
Maximum Weight Capacity	350 lbs.
Maximum Speed	4 mph
Maximum Range	15 miles
Maximum Grade/Incline	6°
Turning Radius	54"
Dimensions	
Overall Length	41.5"
Overall Width	22.5"
Floor to Seat Height	20"-22"
Ground Clearance	3.75"
Power	
Motor	24V x 350W x 3600rpm
Controller	CTE 70A
Brakes	Electromagnetic
Freewheel Mode	Yes
Tires (Flat-free)	
Front	8" x 2"
Rear	8" x 2"
Weights	
Base Weight	77 lbs.
Battery Weight	31 lbs.
Seat Weight	21 lbs.
Total Weight	129 lbs.
Heaviest Piece Weight	46 lbs. (Front Section)

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.