

**Prowler 3310 3-Wheel Heavy Duty Scooter
Order Form/Price List**



Company Name:		Phone:	
Contact Person:		P.O. #	Date:
Bill to:		Ship to:	

Scooter with Seating

- Captain's Seat**
- PROWLER3310MG20CS..... 20" (W) x 20" (D)\$5,311.00
 - PROWLER3310MG22CS..... 22" (W) x 20" (D)\$5,456.50

Batteries & Charger

- 12V x 50AH Batteries (pair) Included
- 5A Battery Charger Included

Accessories

- AB1000 Armrest Bag\$25.00
- AB1120 Backpack\$56.00

- AH1000 Crutch/Cane Holder\$53.00
- AH2000 Cylinder Tank Carrier\$70.00
- AB2100 Rear Basket\$70.00

- AZ3000 Scooter Cover\$72.00
- ST301-COVER 20" Seat Cover\$118.17
- ST306-COVER 22" Seat Cover\$118.17

- AT1000 Scooter Trailer\$204.00

- AZ0060 Cup Holder\$30.00

Specifications

HCPCS	K0899
Performance	
Maximum Weight Capacity	500 lbs.
Maximum Speed	7 mph
Maximum Range	23 miles
Maximum Grade/Incline	7.5°
Turning Radius	50"
Dimensions	
Overall Length	51"
Overall Width	26"
Floor to Seat Height	25"-28"
Ground Clearance	5"
Power	
Motor	24V x 1350W x 4350rpm
Controller	Dynamic Rhino II 120A
Brakes	Electromagnetic & Drum
Freewheel Mode	Yes
Tires (Pneumatic)	
Front	12" x 3.5"
Rear	12" x 3.5"
Weights	
Base Weight	130 lbs.
Battery Weight	66 lbs.
Seat Weight	55 lbs.
Total Weight	251 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.

HCPCS codes are not intended to be, nor should be considered billing or legal advice. Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare Program and should consult an attorney or other advisor to discuss specific situations in further detail.