

Ventura 3-Wheel Midsize/Standard Scooter

Order Form/Price List



Company Name:		Phone:	
Contact Person:	P.O. #	Date:	
Bill to:		Ship to:	

Scooter with Seating

Stadium Seat
 VENTURA318FS Standard, 18" (W) x 17.5" (D).....\$3,112.00

Captain's Seat
 VENTURA318CS Deluxe, 18" (W) x 18" (D)\$3,367.00
 VENTURA320CS Deluxe, 20" (W) x 18" (D)\$3,549.00

Batteries & Charger

12V x 20AH Batteries (pair) Included with STD
 12V x 35AH (U1) Batteries (pair)..... Included with DLX
 4A Battery Charger Included

Accessories

AB1000 Armrest Bag\$25.00
 AB1120 Backpack\$56.00

AH2000 Crutch/Cane Holder\$53.00
 AH2000 Cylinder Tank Carrier\$70.00
 AB2100 Rear Basket\$70.00

AZ3000 Scooter Cover\$72.00
 ST205-COVER 18" Seat Cover\$118.17
 ST301-COVER 20" Seat Cover\$118.17

AT1000 Scooter Trailer\$204.00

AZ0060 Cup Holder\$30.00

Specifications

HCPCS **K0807 (VENTURA318FS Only)**

Performance

Maximum Weight Capacity 350 lbs.
 Maximum Speed 5 mph
 Maximum Range 15 miles (STD) | 22 miles (DLX)
 Maximum Grade/Incline 7.5°
 Turning Radius 53"

Dimensions

Overall Length 48"
 Overall Width 25"
 Floor to Seat Height 19.5"-22.5"
 Ground Clearance 2.5"

Power

Motor 24V x 350W x 5100rpm
 Controller PG S-Drive 70A
 Brakes Electromagnetic
 Freewheel Mode Yes

Tires (Flat-free)

Front 10" x 3"
 Rear 10" x 3"

Weights

Base Weight 88 lbs.
 Battery Weight 28 lbs. (STD) | 46 lbs. (DLX)
 Seat Weight 26 lbs. (STD) | 40 lbs. (DLX)
 Total Weight 142 lbs. (STD) | 174 lbs. (DLX)

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Information contained herein is correct at the time of publication. Drive Medical reserves the right to alter specifications without notice.